
dreambaby®

RAPID RESPONSE CLINICAL THERMOMETER - INSTRUCTIONS
SCHNELLE ANTWORT KLINISCHES THERMOMETER - ANLEITUNG
KLINISCHE THERMOMETER MET KORTE REACTIETIJD - INSTRUCTIES
THERMOMÈTRE CLINIQUE À RÉPONSE RAPIDE - INSTRUCTIONS
TERMOMETRO CLINICO A RISPOSTA RAPIDA - INSTRUZIONI
RESPUESTA RÁPIDA TERMÓMETRO CLÍNICO - INSTRUCCIONES

For the whole family

G338

❶- On/Off Button
❷- LCD Display
❸- Sensor


❶- Ein / Aus-Taste
❷- LCD-Anzeige
❸- Sensor

❶- Aan/Uit Knop
❷- LCD-scherm
❸- Sensor

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❶- Botón encender/apagar
❷- Visualizador LCD
❸- Sensor

READ AND FOLLOW THESE INSTRUCTIONS CAREFULLY AND KEEP FOR FUTURE REFERENCE.

EN Features

- Accurate temperature reading.
- Oral reading in 10 seconds (can vary for underarm and rectal reading).
- Fever Alert - a rapid beep will sound if fever is detected.
- Auto switch off.
- Memory function.
- Flexible tip for comfort.
- Case for hygienic storage.
- Easy to read digital display.
- Contains no mercury.
- Long life battery.
- Celsius and Fahrenheit.
- Temperature range 32.0º – 42.9ºC.
- On / Off button.

To Use

1. Press the On/Off Button next to LCD display.
2. Position thermometer in desired location (mouth, rectum, or armpit).
3. After approximately 10 seconds an alarm will beep and the measured reading will appear on the LCD.

Under the armpit-
Position tip of sensor under the (unclothed) armpit and hold securely till tone is heard. Keep child still as movement may affect accuracy of temperature reading.

Under the tongue-
For children who are two years old or younger, please do not use the devices orally.
Place tip of sensor under the tongue, close mouth and breathe through nose.
Hold securely until beep is heard.
Digital display will show temperature.

Rectally-
Gently insert ONLY tip of sensor into rectum (coat tip with Petroleum Jelly to make insertion easier) and hold securely till tone is heard.

Normal Reading

If the temperature is **below 37.8ºC** a tone will be signalled.

Fever Reading

If the temperature is **above 37.8ºC** a rapid tone will be signalled. If this occurs seek medical advice.

Button

To save power, this thermometer will automatically switch off in 10 minutes.
Turn off the unit by just lightly pressing power.

Memory Function

The next time you use the thermometer the previous temperature reading will show when you turn power on with **Button**.

Changing Celsius and Fahrenheit

1. Ensure device is off.
2. Press and HOLD power **Button**.
 The **Display Panel** will show either "ºF" or "ºC" depending on current setting.
3. Repeat the process to revert back.

To Clean

Clean sensor tip before each use.
Do not submerge the thermometer in cleaning liquid. Use cloth dipped in cold soapy water to wipe clean.
To clean, wipe tip using soft cloth with isopropyl alcohol.
Never dip tip directly into alcohol.
Do not use any other cleaning chemicals.

To Dispose:

Discard this product promptly in accordance with local laws and regulations

To Replace Battery

1. Replace battery when " ▬ " appears in the lower right corner of LCD display.
2. Using a pen, push down on the pogo pin and slide off cover (A).
3. Remove old battery (B). Discard battery lawfully.
4. Replace with new 1.5V DC button type LR41 or SR41, UCC392, or equivalent. Be sure battery is installed with "+" polarity facing up.
5. Slide battery chamber back into place.

Includes: 1x Thermometer, 1x Storage Case

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WARNING EN

CONTAINS BUTTON OR COIN CELL BATTERY. HAZARDOUS IF SWALLOWED. KEEP BATTERIES OUT OF REACH CHILDREN. A BATTERY CAN CAUSE SEVERE OR FATAL INJURIES IN 2 HOURS OR LESS IF IT IS SWALLOWED OR PLACED INSIDE ANY PART OF THE BODY.

- A swallowed Button or Coin Cell Battery can cause internal chemical burns in as little as two hours and lead to death. Due to chemical burns and potential perforation of the esophagus.
- If you think batteries might have been swallowed or placed inside any part of the body, **SEEK IMMEDIATE MEDICAL ATTENTION**.
- If you suspect your child has swallowed or inserted a button battery immediately call the Poisons control hotline and **SEEK IMMEDIATE MEDICAL ATTENTION**.
- Examine devices and make sure the battery compartment is correctly secured. Do not use if compartment is not secure.
- Tell others about the risk associated with button batteries and how to keep their children safe.
- Do not dispose of batteries in fire, or expose to extreme heat, it may explode.

Caution: To prevent battery leaks, which can burn skin and eyes:
•Remove battery when storing product for a long time.
•Discard old batteries promptly in accordance with local laws and regulations.

Precaution: This thermometer is a precise electronic device, designed to provide accurate temperature readings.
To remain reliable:
- NEVER immerse into hot or boiling water or other liquids.
- NEVER chew or bite any part of the thermometer.
- NEVER bend or subject the thermometer to any sudden, heavy shocks.
- NEVER subject thermometer to direct sunlight, extreme temperatures, humidity, dust etc.
- NEVER dismantle any part of the thermometer, apart from the battery cover, when replacing the battery.
- Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- Note: Use of the probe cover may result in a 0.1°C discrepancy from actual temperature.
- Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- If the unit has been stored at temperatures over 5°C~40°C, leave it in 5°C~40°C ambient temperature for about 15 minutes before using it.

This item is not a toy. Keep out of reach of children when not in use. Children should be under close supervision at all times while using this product.

Certification: This thermometer is designed and manufactured complying with relevant requirements. It has passed EMC emission test, immunity test according to EN60601-1-2 and biocompatibility test according to Medical Device Directive annex 1. It is rated BF for protection against electric shock, in accordance with EN60601-1. This appliance conforms to the following standards: ISO 80601-2-56 Medical electrical equipment —Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment — Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2 (EMC), IEC/EN60601-1 (Safety) standards. Manufacturer is ISO 13485 certified.

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©2026 Tee-Zed Products Ltd
Tee-Zed Products Pty Ltd
Tee-Zed Products Ltd
Tee-Zed Products (EU) Ltd