dreambaby®

INSTRUCTIONS • INSTRUCCIONES

Model • Modelo: L341 (DET-306)

Infrared Forehead Thermometer

2x AAA batteries (included)

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2x AAA pilas (incluidas)

Introduction

The Infrared Forehead Thermometer is a device capable of measuring people's body temperature by detecting the intensity of infrared light emitted from the forehead.
It converts the measured heat into a temperature reading displayed on the LCD.
The infrared forehead thermometer is intended for the intermittent measurement of human body temperature from the skin surface of forehead by people of all ages.
When properly used, it will quickly assess your temperature in an accurate manner.

Please read all instructions carefully and thoroughly before using this product.

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
WARNING:
- There is no gender or age limitation for using infrared forehead thermometer.
- This thermometer is intended for home use only.
- Use of this Forehead thermometer is not intended as a substitute for consultation with your physician.
- Do not allow children to take their temperatures unsupervised, some parts are small enough to be swallowed.
- Never immerse this device in water or other liquids (not waterproof).
- Do not modify this equipment without authorization of manufacturer.
- Do not expose the thermometer to temperature extremes (below -25ºC / -13ºF or over 55ºC / 131ºF nor excessive humidity (>95%RH).
- Keep the battery away from children.
- Remove battery from the device when not in operation for a long time.

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
Forehead Thermometer Advantages

Infrared Forehead Thermometer measures core body temperature, which is the temperature of a body's vital organs (See Figure 1).
This thermometer is designed to measure the temperature of the skin surface over the temporal artery, a major artery of the head.
The temporal artery is connected to the heart via the carotid artery, directly leading from the aorta, the main trunk of the arterial system.
It offers constant blood flow. Therefore, body temperature changes are reflected sooner in the forehead than they are in other parts of the body such as oral, rectal and underarm.

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FIG. 1

Temperature Taking Hints

To ensure that the reading always reflects the body temperature accurately, you need to take account of the following factors which may affect an accurate reading.
1. It is important to know each individual's normal temperature when they are well. This is the only way to accurately diagnose a fever. To determine normal temperature, take multiple readings when healthy. Re-measure with a standard digital thermometer for confirmation.
2. Users must be inside for 30 minutes before taking a measurement.
 Note: Users and the thermometer should be in the same ambient temperature for at least 10 minutes before taking a reading.
3. Users should not drink, eat, or be physically active such as bathing, showering, shampooing and hair drying before/while taking the measurement. Remove hat and wait 10 minutes before taking a reading.
4. Oils or cosmetics on the forehead may give a lower temperature reading than the actual one. Remove dirt from the forehead before taking a measurement. Wait at least 10 minutes after washing the forehead area before taking a reading.
5. Holding a hand on the forehead for any length of time will affect the temperature reading.
6. Do not take temperature over scar tissue, open sores or abrasions.
7. Do not use the thermometer on a perspiring or sweating forehead, as this may affect the reading.
8. Don't take a measurement while or immediately after nursing a baby.
9. This thermometer is intended for home and clinical use.
10. Do not take temperatures with this thermometer near places that are very hot, such as fireplaces and stoves.
11. The probe window of the thermometer is the most delicate part of the device. Do not touch the probe window. The accuracy of the reading may be affected if the probe window is damaged or dirty.
12. If the thermometer is stored in a significantly different environment than testing location, place it in the testing location for approximately 30 minutes prior to use.
13. It is not intended for use in the oxygen rich environment and presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.

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
Clock Setting

When using thermometer for the first time, please set the parameters of the thermometer.
With the thermometer off, press and hold ⏻ BUTTON to enter into setting mode for 3 seconds.

1. Set the time format
The device can display the time in either an AM/PM (12-hour) or a 24:00 (24-hour) format. Press and release START BUTTON to select the format. With the preferred time format on the display, press ⏻ BUTTON, the Hour figure is flashing automatically.

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2. Set the hour
Press and release the START BUTTON to advance one hour until the correct hour appears. After the hour is set, press ⏻ BUTTON, the Minute figure is flashing automatically.

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3. Set the minute
Press and release the START BUTTON to advance one minute until the correct minute appears. After the minute is set, press ⏻ BUTTON, the Year figure is flashing automatically.

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4. Set the year
Press and release the START BUTTON to advance one year until the correct year appears. After the year is set, press ⏻ BUTTON, the Month figure will appear.

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5. Set the month
Press and release the START BUTTON to advance one month until the correct month appears. After the month is set, press ⏻ BUTTON, the Date figure is flashing automatically.

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6. Set the date
Press and release the START BUTTON to advance one day until the correct month appears. After the day is set, press ⏻ BUTTON to exit the setting mode.

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To Use

▶ To measure forehead temperature:
1. Press the ⏻ BUTTON, the display is activated to show all segments.
2. After self-checking Figure 2 appears on the display screen with beeps, so you can start a new measurement. Aim the thermometer at the centre of the forehead with a distance of 2 to 3cm (See figure 3) and then press the START BUTTON.
 Note: Do not remove the thermometer from the forehead before hearing the final beep. A waiting period of 3 seconds between testing is recommended to ensure complete accuracy.
3. Read the temperature on the display.
4. Press the ⏻ BUTTON to turn off.

FIG. 2

FIG. 3

▶ Alert - Backlight:
In Forehead mode:
1. The display will be lighted GREEN for 3 seconds with a happy face ☺ when the unit is ready for measurement and a measurement is completed with a reading less than 37.3ºC (99.1ºF).
2. The display will be lighted YELLOW for 3 seconds with a happy face ☺ when a measurement is completed with a reading less than 37.8ºC (100.0ºF).
3. The display will be lighted RED for 3 seconds with a sad face ☹ when a measurement is completed with a reading equal to or higher than 37.8ºC (100.0ºF).

▶ After measurement:
1. Power off: Device will automatically shut off if left idle for more than 1 minute to extend battery life.
2. Clean the probe after each use to ensure an accurate reading and avoid cross contamination. (See the section of care and cleaning for details.)

▶ How to change forehead mode and object mode:
1. You can press 👤/🏠 BUTTON to switch the mode between forehead mode and object mode.

▶ To measure object temperature:
1. Press the ⏻ BUTTON to turn on the thermometer, you can take the object temperature after hearing two beep sounds (see figure 4).
2. Aim the thermometer at the centre of the object you want to measure with a distance of 1 to 2cm.
3. Press the START BUTTON and then read the temperature on the display.
4. Press the ⏻ BUTTON to turn off.


FIG. 4

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Memory Mode

1. The Memory Mode can be accessed ether in forehead mode or object mode: When the thermometer has been turned on and followed by Figure 2/4 or finished testing, press and hold the ⏻ BUTTON for three seconds. The letter M will appear in the centre right corner of the display. (See Figure 5)
2. The thermometer will automatically memorize the last 10 temperature readings. Each memory also records the measurement date/ time/mode icon. Each time the ⏻ BUTTON is pressed, the screen displays past readings that correspond with a number 1-10. The number 1 reflects the most recent reading, while the number 10 reveals the oldest reading stored in memory. (See Figure 6)
3. In the memory mode, 👤mark or 🏠 mark always exist. The user can press the START BUTTON to take new measurements.

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
FIG. 7

Selecting A Temperature Scale

1. Temperature readings are available in the Celsius (ºC) or Fahrenheit (ºF) scale.
2. With the thermometer off, press and hold the START BUTTON for 3 seconds to enter into unit changing mode.
3. Press and release START BUTTON to select the unit.
4. When the preferred unit on the display, press ⏻ BUTTON to exit the unit changing mode. (See Figure 7)

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Care and Cleaning

1. The probe window must be kept clean, dry, and undamaged at all times to ensure accurate readings. The accuracy of temperature readings can be affected by damage to the probe window, or the presence of dirt, fingerprints, earwax, dust and other soiling compounds on the probe window. Degraded sensors can degrade performance or cause other problems.
2. The probe window is the most delicate part of the thermometer. Use a soft cloth slightly moistened with a 75% isopropyl alcohol solution to disinfect probe window and the thermometer. Do not use abrasive cleaners. After cleaning, allow at least 10 minutes drying time before taking temperatures. **Note:** Do not use any chemical other than isopropyl alcohol to clean the probe window.
3. Use a soft, dry cloth to clean the thermometer display and exterior.
4. The thermometer is not waterproof. Do not submerge the unit in water when cleaning.
5. Store the thermometer in a dry location, free from dust and contamination and away from direct sunlight.
6. Strong electromagnetic fields may interfere with the proper operation of the thermometer. The device needs special pre-cautions regarding EMC according to the EMC information.
7. Put the thermometer back to the original packaging after using.

Battery Replacement

1. Replace battery when ⌧ appears in the upper right corner of LCD display (see Figure 8).
2. Slide battery cover down as shown in Figure 9.
3. Remove battery and install 2 new AAA alkaline (or 2 x AAA 1.5V rechargeable) batteries as shown in Figure 10.
4. Slide battery cover back on.

FIG. 8

FIG. 9

FIG. 10

- Ensure the batteries are installed correctly according to polarity (+ and -).
- Do not mix old and new batteries, different brands or types of batteries, such as alkaline, carbon-zinc, or rechargeable batteries.
- Remove and immediately recycle or dispose of batteries from equipment not used for an extended period of time according to local regulations.
- Always completely secure the battery compartment. If the battery compartment does not close securely, stop using the product, remove the batteries, and keep them away from children.

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Electromagnetic Compatibility Information

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

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Warning:
1. The Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.)
2. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the [ME EQUIPMENT or ME SYSTEM], including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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For transmitters rated at a maximum output power not listed above, the recommended separation distance "**d**" in metres (**m**) can be estimated using the equation applicable to the frequency of the transmitter, where "**P**" is the maximum output power rating of the transmitter in watts (**W**) according to the transmitter manufacturer.

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NOTE 1 - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

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NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Certification

This appliance conforms to the following standards:

ISO 80601-2-56 Medical electrical equipment

- Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, IEC 60601-1-11 Medical electrical equipment
- Part 1-11: General requirements for basic safety and essential performance
- Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of IEC 60601-1-2 (EMC), IEC 60601-1(Safety) standards.

And the manufacturer is ISO 13485 certified.

Calibration

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the use instruction, periodic readjustment is not required. However, we recommend checking calibration every two year or whenever clinical accuracy of the thermometer is in question. Please send the complete device to the dealer or manufacturer.

The above recommendations do not supersede the legal requirements. The user must always comply with legal requirement for the control of the measurement, functionality, and accuracy of the device which are required by the scope of relevant laws, directives or ordinances where the device is used.

A clinical summary and procedures for checking calibration are available upon request. (Turn of the thermometer and press the power button for a long time until entering into calibrate mode, software version will be displayed).

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TEE-ZED PRODUCTS, LLC ACCEPTS NO RESPONSIBILITY FOR ANY LOSS OR DAMAGE SUFFERED BY ANY PERSON AS A RESULT OF THE USE OR MISUSE OF THESE GOODS. IN THE CASE OF ANY DEFECT IN THE GOODS, TEE-ZED'S LIABILITY SHALL BE LIMITED SOLELY TO REPLACEMENT OF THE DEFECTIVE GOODS.

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Tee-Zed Products Inc.
401 Norfinch Drive,
North York Ontario, M3N 1Y7
CA Toll Free/Sans Frais
1-877-957-2229
info@dreambaby.ca
www.dreambaby.ca

AUSTRALIA
info@tee-zed.com.au
www.dreambaby.com.au

#dreambabysafety
YouTube
dreambabytv

Made in Shenzhen, China
Hecho en Shenzhen, China
©2026 Tee-Zed Products LLC
Manufactured by • Fabricado por: JOYTECH HEALTHCARE CO., LTD.
No.365, Wuzhou Road, Yuhang Economic Development Zone,
311100, Hangzhou, China

X.04

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Tee-Zed Products, LLC